



# Digital Stethoscopes in the Age of AI and Connected Medicine: Clinical Implications for Practice

Andrès E<sup>1,2</sup>, El Hassani HA<sup>3</sup>, Talha S<sup>2,4</sup>, Lavigne T<sup>2,5</sup>, Terrade JE<sup>1</sup>, Jannot X<sup>1</sup> and Lorenzo-Villalba N<sup>1</sup>

<sup>1</sup>Department of Internal Medicine, Strasbourg University Hospitals (HUS), France

<sup>2</sup>Mitochondria, Oxidative Stress and Muscle Protection Research Unit (MSP), Faculty of Medicine, University of Strasbourg, France

<sup>3</sup>Laboratory of Nanomedicine, Imaging and Therapeutics, University of Technology of Belfort-Montbéliard (UTBM), France

<sup>4</sup>Physiology and Functional Exploration Unit, Strasbourg University Hospitals (HUS), France

<sup>5</sup>Public Health Services, Strasbourg University Hospitals (HUS), France

## Abstract

For over two centuries, the stethoscope has been a cornerstone of clinical medicine, yet traditional acoustic models are limited by environmental noise, user experience, and the subjective nature of auscultation. Digital stethoscopes, augmented by signal processing, visualization tools, and Artificial Intelligence (AI), represent a transformative evolution in auscultatory practice. These devices convert acoustic signals into digital data, enabling amplification, noise reduction, and objective analysis through phonograms and spectrograms. AI integration enhances diagnostic accuracy by detecting and classifying cardiac and respiratory abnormalities, while connectivity features facilitate telemedicine, remote monitoring, and longitudinal patient care. Clinical studies consistently demonstrate the superiority of digital over acoustic stethoscopes in detecting subtle heart murmurs, lung sounds, and arrhythmias, particularly in complex or noisy settings. Beyond clinical diagnostics, digital stethoscopes serve as educational tools, providing visualized, recordable, and analyzable auscultatory data. As part of the emerging Medicine 3.0 ecosystem, these devices align with predictive, preventive, personalized, participatory, and precision (5P) medicine, offering new opportunities for global health equity, chronic disease management, and medical education. This review synthesizes current evidence on the technology, clinical performance, AI applications, and telemedicine integration of digital stethoscopes, highlighting their potential to redefine auscultation in the 21st century.

**Keywords:** Digital stethoscope; Auscultation; Artificial intelligence; Telemedicine; Phonogram; Spectrogram; Cardiac diagnosis; Respiratory diagnosis; Medicine 3.0; Precision medicine.

## Key Points

- **Digital Transformation of Auscultation:** Digital stethoscopes convert acoustic signals into digital data, enabling amplification, noise reduction, and objective visualization through phonograms and spectrograms.
- **AI-Enhanced Diagnostics:** Integrated artificial intelligence improves the detection and classification of cardiac and respiratory abnormalities, enhancing diagnostic accuracy and reducing inter-observer variability.
- **Telemedicine Integration:** Connectivity features allow remote auscultation, expert consultation, and longitudinal patient monitoring, supporting healthcare delivery in underserved and remote regions.
- **Clinical Superiority:** Evidence from comparative studies demonstrates higher sensitivity and specificity for subtle heart murmurs, arrhythmias, and lung sounds compared to traditional acoustic stethoscopes, especially in noisy or complex clinical settings.
- **Educational and Research Utility:** Digital stethoscopes provide recordable, visualized, and analyzable auscultatory data, enhancing medical education, training, and research in cardiopulmonary medicine.
- **Alignment with Medicine 3.0:** By enabling predictive, preventive, personalized, participatory, and precision (5P) medicine, digital stethoscopes contribute to improved patient outcomes, chronic disease management, and global health equity.

## INTRODUCTION

For over two centuries, the stethoscope has symbolized the essence of clinical medicine—a bridge between physician and patient, ear and heart, observation and intuition [1]. Introduced by René Laennec in 1816, the acoustic stethoscope was not merely a diagnostic tool but a transformation of clinical method, allowing physicians to “see” within the

body through sound. Despite its iconic status, the acoustic stethoscope has remained largely unchanged while medicine has progressed into the digital age. Today, digital stethoscopes stand at the forefront of a technological revolution in auscultation [2]. By converting acoustic signals into digital data, providing visualization via phonograms and spectrograms, integrating Artificial Intelligence (AI) for automated analysis, and enabling telemedicine applications, these devices enhance diagnostic accuracy, reduce human error, and extend access to care across diverse settings. Beyond technical advancement, this evolution invites a rethinking of what it means to “listen” in modern medicine. This review aims to provide a comprehensive overview of digital stethoscopes, including their technological components, signal-processing capabilities, clinical performance compared with acoustic models, AI integration, telemedicine applications, and future directions in the era of connected, data-driven healthcare.

## TERMINOLOGY AND DEVICE CLASSIFICATION

The terminology surrounding stethoscopes has evolved to reflect the expanding technological diversity of modern auscultatory devices [3].

### Acoustic Stethoscopes

Remain the traditional standard. They consist of a chest piece—

**Submitted:** 19 November 2025 | **Accepted:** 06 December 2025 | **Published:** 08 December 2025

**\*Corresponding author:** Emmanuel ANDRES, Head of the Department of Internal Medicine, Haute-pierre Hospital, Strasbourg University Hospitals 1, Avenue Molière, 67091 STRASBOURG Cedex

**Copyright:** © 2025 Andrès E, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

**Citation:** Andrès E, El Hassani HA, Talha S, Lavigne T, Terrade JE et al. (2025) Digital Stethoscopes in the Age of AI and Connected Medicine: Clinical Implications for Practice. SM J Pulm Med 7: 9.



typically with a diaphragm and/or bell—connected to flexible tubing transmitting sound to earpieces. Operating purely on mechanical principles, they provide reliable auscultation without electronic amplification or signal processing, making them simple, robust, and widely accessible.

## Electronic Stethoscopes

Introduce an initial layer of technology. A built-in microphone converts mechanical sound waves into electronic signals, which are amplified and filtered to improve audibility, particularly in noisy environments. While the output is still analog, internal digitization enhances clarity, reduces distortion, and facilitates sound storage. Additional features such as volume control and selective noise filtering make electronic stethoscopes particularly suitable for busy clinical settings.

## Digital Stethoscopes

Represent a significant advance by converting sound into fully digital data via analog-to-digital converters (ADCs) [3]. This enables visualization of sounds through phonograms and spectrograms, long-term storage, advanced filtering, and wireless transmission for telemedicine. Digital models frequently incorporate noise reduction algorithms to enhance subtle cardiac or pulmonary sounds, supporting more precise detection of abnormalities.

## Smart or AI-Enabled Stethoscopes

Integrate artificial intelligence algorithms, either embedded or cloud-based, to assist clinicians in detecting and classifying pathological sounds, such as murmurs, wheezes, and crackles [3]. AI enhances diagnostic precision, reduces inter-observer variability, and can serve as a training tool for clinicians by providing feedback on auscultation technique.

## Numerical Stethoscopes

Quantify auscultatory data, converting sounds into measurable metrics such as frequency, amplitude, and heart or respiratory rates [3,4]. When combined with AI, these devices offer objective, data-driven support for diagnosis, longitudinal monitoring, and chronic disease management. Although these categories overlap, distinguishing between them is essential for evaluating regulatory requirements, clinical performance, and application in patient care. As innovation continues, the boundaries between electronic, digital, and AI-enabled stethoscopes are increasingly blurred, redefining the practice of auscultation in modern medicine.

## TECHNOLOGICAL OVERVIEW: FROM SOUND TO SIGNAL

Digital stethoscopes combine several advanced technological components to augment traditional auscultation, overcoming limitations inherent to acoustic devices [3,4]. Understanding these elements is essential to appreciating their clinical utility and diagnostic potential.

### Transducers and Microphones

The transducer is the initial interface with the patient, housing microphones designed to capture a broader range of frequencies than conventional diaphragms. Some models integrate dual microphones optimized for low-frequency heart sounds (bell) and high-frequency lung sounds (diaphragm), electronically merging signals for comprehensive auscultation. Modern microphones minimize ambient noise while enhancing subtle physiological sounds.

### Amplification and Filtering

Unlike acoustic stethoscopes, where sound intensity diminishes along the tubing, digital stethoscopes use electronic amplification to enhance volume. Clinicians can adjust sound levels to detect faint murmurs or subtle respiratory sounds. Digital signal processing (DSP) further applies frequency-specific filters—low, mid, and high—to isolate

clinically relevant sounds, improving clarity in noisy environments or in patients with thick chest walls.

### Analog-to-Digital Conversion (ADC)

Captured sounds are digitized using ADCs, typically with 16–24 bit resolution. The sampling rate and bit depth determine the fidelity of the digital signal, preserving subtle variations in heart and lung sounds [3–5]. High-resolution digitization ensures precise representation, enabling advanced analysis and storage.

### Digital Signal Processing (DSP)

DSP algorithms enhance the digitized audio in real time, applying filtering, compression, and noise reduction. DSP also optimizes signals for storage, transmission, and visualization. Spectral analysis generates phonograms and spectrograms, providing graphical representations of sound amplitude and frequency over time [6,7]. These visualizations facilitate detection of abnormal patterns and support objective clinical assessment.

### Data Storage and Connectivity

Many digital stethoscopes include onboard memory for storing auscultation recordings. Wireless connectivity via Bluetooth or Wi-Fi allows transmission to Electronic Health Records (EHRs), remote devices, or telemedicine platforms, enabling expert consultation without physical presence. Real-time streaming supports synchronous tele-auscultation, improving accessibility in remote or underserved areas.

### Artificial Intelligence Integration

Advanced digital stethoscopes incorporate AI modules—either embedded or cloud-based—to assist in pattern recognition and diagnostic suggestions. AI algorithms can detect heart murmurs, wheezes, crackles, arrhythmias, and other pathological sounds, enhancing clinician accuracy, reducing inter-observer variability, and supporting training [4].

### Power Source

Most devices rely on rechargeable lithium-ion batteries, offering several hours of continuous use. Battery life is a critical consideration to ensure reliable performance in clinical practice.

### Comparative Advantages


Table 1 summarizes the differences between acoustic and digital stethoscopes, highlighting enhancements in sound amplification, frequency range, noise reduction, data storage, and connectivity [3,4]. These features collectively represent a transformative step forward in auscultatory medicine, enabling clinicians to combine auditory, visual, and data-driven approaches for improved patient care.

## THE INFORMATIVE POWER OF SPECTROGRAMS AND PHONOGRAMS: LISTENING WITH THE EYES

The digitization of auscultatory sounds has enabled the creation of visual representations, providing objective, quantifiable data that complement traditional auditory assessment [4] (Figure 1). This innovation allows clinicians to “see” sounds, enhancing diagnostic precision, documentation, and education.

### Phonograms

Phonograms graphically display sound amplitude over time, with the x-axis representing time and the y-axis showing sound intensity. These representations allow detailed analysis of timing, duration, and intensity of heart and respiratory sounds [7]. For instance, the crescendo-decrescendo pattern of an aortic stenosis murmur or the flat holosystolic profile of mitral regurgitation can be clearly visualized. Phonograms also



**Figure 1:** Educational visualization of cardiac sounds as a cardio-phonogram and cardio-spectrogram.

Legend: Left: image generated by artificial intelligence. Right: image from auscultation of a patient with aortic stenosis - the latter is indicated by a white arrow during a cardiac cycle materialized by B1-B2 (personal data; ASAP project).

**Table 1:** Comparison between Acoustic and Digital stethoscopes.

Feature	Acoustic Stethoscope	Digital Stethoscope
Sound Transmission	Relies on mechanical sound transmission through tubing	Uses digital sensors to capture and amplify sound
Sound Quality	Prone to background noise interference and attenuation	Enhanced sound clarity with noise reduction features
Frequency Range	Typically limited to human audible range (20-1,000 Hz)	Wider frequency range, often covering 20 Hz to 2,000 Hz+
Amplification	No amplification, relies on the user's ear sensitivity	Digital amplification allows louder and clearer sounds
Battery Life	No batteries required	Requires batteries, with varying life depending on usage
Recording Capability	No recording capability	Can record and store auscultation sounds for later analysis
Connectivity	No digital connectivity	Often integrates with smartphones, computers, or telemedicine platforms
Data Storage/Analysis	None	Can analyze and store data, sometimes with advanced algorithms or AI
Size/Weight	Generally lightweight and portable	May be slightly heavier due to digital components
Price	Typically more affordable	Higher cost due to advanced technology and features
User Experience	Simple, familiar design	Requires some technical knowledge to operate effectively
Use in Telemedicine	Limited use without additional equipment	Ideal for integration into telemedicine and remote diagnostics
Maintenance	Low maintenance, just tubing and earpiece care	Requires software updates and occasional battery replacement
Durability	Long-lasting with proper care	May require more frequent servicing due to digital components
Customization	Limited customization options	Can be customized with different sound filters and settings



facilitate differentiation of normal heart sounds (S1, S2) from abnormal sounds such as S3, S4, clicks, and snaps. This visual documentation improves communication among healthcare providers and enhances clinical decision-making (Figure 1).

### Spectrograms

Spectrograms provide frequency-based visualizations, with the x-axis as time, the y-axis as frequency (Hz), and color intensity representing amplitude. This allows clinicians to analyze harmonic content and energy distribution across frequencies, which is invaluable in distinguishing sounds with similar temporal patterns but differing frequency profiles. For example, spectrograms can help differentiate wheezes from crackles or classify various cardiac murmurs by their unique frequency signatures [7] (Figure 1).

### Clinical Utility

Visual representations enable longitudinal monitoring of disease progression, response to treatment, and objective evaluation of interventions [4]. They also serve as powerful educational tools, allowing medical students and trainees to learn auscultation with visual feedback, improving pattern recognition and diagnostic accuracy.

### Telemedicine Integration

Phonograms and spectrograms can be transmitted remotely for expert review, facilitating tele-auscultation and collaborative diagnosis. This is particularly valuable in resource-limited settings or for patients in rural areas where access to specialists is limited.

### Evidence

Recent studies indicate that incorporating visual representations of heart and lung sounds improves inter-observer agreement in murmur classification and supports remote education and training [7]. By combining auditory and visual information, clinicians gain a more precise and reproducible understanding of auscultatory findings. In summary, phonograms and spectrograms transform auscultation from a purely auditory skill into a multimodal, data-driven process, enhancing diagnostic accuracy, clinical education, and remote healthcare capabilities.

### THE AUGMENTING ROLE OF ARTIFICIAL INTELLIGENCE IN AUSCULTATION

The integration of Artificial Intelligence (AI) into digital stethoscopes is transforming auscultation by enabling automated analysis of cardiac

**Table 2:** Overview of Commercial Digital Stethoscopes in USA and Europe

Brand/Model	Country Available	Features	Sound Quality	Battery Life	Connectivity	Price Range (USD)	Notable Clinical Features
3M Littmann 3200	USA, Europe	Bluetooth, noise reduction, audio recording	Excellent, clear sound, enhanced with ambient noise reduction	Up to 60 hours	Bluetooth, connects to mobile app	\$300 - \$400	Advanced filtering, integrated app for sound visualization, voice recording
Eko CORE Digital Stethoscope	USA, Europe	Bluetooth, amplifies sounds, app integration	Excellent, customizable filters	Up to 8 hours	Bluetooth, Eko software app	\$250 - \$350	Real-time heart and lung sound visualization, telemedicine integration
Omron Stethoscope (Digital)	USA, Europe	Noise reduction, digital sound amplification	Good, enhanced with noise cancellation	10 - 15 hours	None	\$100 - \$150	Simplified design for daily use, battery-operated, less customization
iHealth Stethoscope	USA, Europe	Amplification, mobile app integration	Clear, moderate amplification	12 - 16 hours	Bluetooth, iHealth app	\$150 - \$200	Audio capture for later analysis, compatible with iOS/Android apps
Thinklabs One	USA, Europe	High amplification, small and portable design	Excellent, very high amplification	Up to 25 hours	None	\$500 - \$600	Exceptional amplification, ideal for difficult-to-hear patients, compact design
Welch Allyn Harvey DLX	USA, Europe	Digital frequency tuning, amplification	High, customizable sound profiles	20 - 25 hours	None	\$250 - \$350	High-quality acoustic and digital sound, multi-frequency adjustments
Amplified Stethoscope by ADC	USA, Europe	Sound amplification, flexible tubing	Moderate, amplified for clarity	12 - 15 hours	None	\$100 - \$175	Simple amplification, affordable price, useful for general auscultation
Littmann Electronic Stethoscope Model 3100	USA, Europe	Digital amplification, noise reduction, dual head design	High-quality with clear sound	Up to 20 hours	Bluetooth, integrates with mobile app	\$250 - \$350	Acoustic and digital dual head design, customizable sound filters, digital recording



and pulmonary sounds, enhancing both accuracy and efficiency in clinical practice [3,4]. AI acts as an augmentative “second ear,” supporting clinicians in detecting subtle or complex abnormalities that may be challenging to identify with conventional auscultation.

### Automated Detection of Cardiac Abnormalities

AI algorithms trained on extensive datasets of normal and pathological heart sounds can detect murmurs, arrhythmias, atrial fibrillation, and valvular heart diseases with high precision [3,4]. These algorithms analyze timing, amplitude, and spectral patterns, providing real-time alerts for subtle anomalies that may otherwise go unnoticed, particularly in early-stage disease or noisy clinical environments.

### Automated Analysis of Respiratory Sounds

AI also enhances the identification of abnormal lung sounds, including wheezes, crackles, and pleural rubs, associated with conditions such as pneumonia, COPD, asthma, and pleural effusions [3,4]. Automated analysis supports early detection, continuous monitoring, and timely intervention, improving patient outcomes.

**Risk Stratification and Prognosis:** Advanced AI models can predict the risk of future cardiac events by analyzing nuanced acoustic features. This predictive capability enables early intervention and personalized treatment planning, shifting auscultation from reactive diagnosis to proactive risk management [3,4].

### Reducing Inter-Observer Variability

By providing objective, data-driven analysis, AI reduces inter-observer variability in auscultation, particularly benefiting less experienced clinicians. Studies have shown that AI-assisted stethoscopes improve diagnostic confidence and consistency, facilitating standardized clinical assessments across diverse healthcare settings [3,4].

### Big Data Integration

Digital auscultation generates large volumes of structured and semi-structured data, including phonocardiograms, spectrograms, and respiratory cycles. AI can aggregate and analyze these data across populations, identifying longitudinal patterns and subtle abnormalities that may be imperceptible to the human ear [3,4,7]. Machine learning models, including convolutional and recurrent neural networks, have demonstrated remarkable accuracy in classifying heart murmurs, detecting atrial fibrillation, and distinguishing between different lung sounds.

### Evidence

A 2023 multicenter study analyzed over 1.2 million auscultation recordings and found that a deep learning model outperformed general practitioners in detecting pathological murmurs (AUC 0.93 vs. 0.79), highlighting the potential of AI to enhance diagnostic precision and clinical decision-making [8].

### Clinical Considerations

While AI significantly augments auscultation, it should complement, not replace, clinical judgment. Results must be interpreted within the broader clinical context, integrating patient history, physical examination, and additional diagnostic tests. In conclusion, AI-powered digital stethoscopes represent a paradigm shift in clinical auscultation, offering improved detection, risk stratification, and standardization, while supporting telemedicine, longitudinal monitoring, and evidence-based decision-making.

## MARKET OVERVIEW AND MODELS OF DIGITAL STETHOSCOPES

The market for digital stethoscopes has expanded rapidly, driven by advances in sound amplification, signal processing, connectivity, and AI integration (4). Modern devices are designed to meet diverse clinical needs, from routine auscultation to specialized cardiology and pulmonology assessments, and increasingly support telemedicine applications.

### Key Features and Functionalities

Digital stethoscopes differentiate themselves from traditional acoustic models through several capabilities:

**Sound Amplification:** Enhances low-intensity cardiac and pulmonary sounds, improving detection of subtle murmurs or faint respiratory noises.

**Noise Reduction:** Active filtering isolates clinically relevant sounds from ambient noise, particularly valuable in emergency departments, ICUs, or outpatient settings.

**Connectivity:** Wireless options such as Bluetooth and Wi-Fi enable seamless integration with mobile devices, telemedicine platforms, and electronic health records (EHRs).

**Data Storage and Sharing:** Onboard memory allows storage of auscultation recordings for longitudinal monitoring, review, or expert consultation.

**AI-Driven Diagnostics:** Certain models include embedded or cloud-based AI algorithms for real-time detection and classification of heart and lung sounds, including murmurs, arrhythmias, wheezes, and crackles.

### Representative Models in the USA and Europe (including France)

A variety of commercial devices now provide combinations of the above features:

**Littmann 3200:** Wireless transmission, sound amplification, and digital storage, widely used in cardiology and general practice.

**Thinklabs One:** Compact design with high-fidelity amplification and recording capabilities, suitable for pediatrics.

**Eko Core:** Integrates AI-assisted murmur detection and ECG functionality, bridging auscultation with cardiac monitoring.

**Stethoscope:** Consumer-friendly design with Bluetooth connectivity for telemedicine, supporting remote auscultation and patient engagement.

### Clinical Applications

These devices have demonstrated utility in improving auscultation in noisy environments, enhancing diagnostic accuracy, supporting teleconsultations, and providing training opportunities for medical students and residents. As connectivity and AI capabilities expand, digital stethoscopes are increasingly becoming essential tools for both clinical practice and telehealth workflows [4].

### Considerations for Clinicians

Selection of a digital stethoscope should be guided by clinical needs, including patient population (adult, pediatric), frequency of remote consultations, AI diagnostic support, battery life, and device ergonomics. While devices differ in features, all provide enhanced capabilities compared to acoustic stethoscopes, particularly in precision, documentation, and telemedicine integration.

The evolving market reflects the growing role of digital stethoscopes as central components of modern, connected healthcare. By combining traditional auscultation with amplification, visualization, AI support, and telemedicine integration, these devices offer a comprehensive approach



to cardiac and pulmonary assessment, redefining the clinical experience of listening to the human body (Table 3 and Table 4).

## RESULTS OF STUDIES COMPARING ACOUSTIC AND DIGITAL STETHOSCOPIES

A growing body of evidence has evaluated the comparative performance of digital versus traditional acoustic stethoscopes, highlighting the clinical advantages of digital devices in diverse settings.

### Superior Detection of Low-Intensity Sounds

Digital stethoscopes with amplification capabilities significantly improve the detection of faint heart murmurs (grade 1–2) and subtle respiratory sounds [3,4]. This benefit is particularly evident in noisy environments or in patients with thick chest walls, where acoustic

stethoscopes may fail to capture low-intensity signals.

### Enhanced Diagnostic Accuracy through Visualization

Devices equipped with phonogram and spectrogram visualization allow clinicians to objectively characterize heart and lung sounds [7]. Visual representations assist in distinguishing S1–S2 from S3–S4, murmurs, clicks, and abnormal respiratory sounds, enhancing diagnostic precision, especially for less experienced practitioners.

### AI-Assisted Screening and Diagnosis

AI algorithms integrated into digital stethoscopes have shown high sensitivity and specificity in detecting conditions such as atrial fibrillation, heart murmurs, and pulmonary abnormalities [3,8]. These algorithms act as a “second ear,” alerting clinicians to subtle abnormalities that may be

**Table 3:** Digital stethoscopes that integrate Artificial Intelligence (AI) for enhanced clinical use

Brand/Model	AI Features	Sound Quality	Connectivity	Key Clinical Features	Price Range (USD)
Eko CORE Digital Stethoscope	AI-driven algorithms for heart and lung sound analysis	Excellent, customizable filters	Bluetooth, integrates with Eko software app	Real-time visualization of heart and lung sounds, AI-based arrhythmia detection, telemedicine integration	\$250 - \$350
3M Littmann 3200 - CORE	AI-assisted noise reduction and sound filtering for clearer auscultation	Excellent, clear with noise reduction	Bluetooth, mobile app integration	AI-driven noise reduction, custom sound profiles, sound visualization	\$300 - \$400
Eko DUO	Combines digital amplification with AI-based ECG analysis	Excellent, enhanced with AI filtering	Bluetooth, integrates with Eko software	Dual mode (stethoscope + ECG), AI for detecting arrhythmias, heart murmurs	\$450 - \$600
Thinklabs One with Thinklabs Digital	AI-assisted sound analysis for enhanced diagnostics	Excellent, very high amplification	None (does not directly integrate with apps)	AI-driven noise reduction, enhanced amplification for tough-to-hear patients	\$500 - \$600
Steth IO	AI-powered sound analysis with real-time diagnostics	Clear, amplified sounds	Bluetooth, mobile app integration	AI-assisted lung and heart sound analysis, remote diagnosis capabilities	\$300 - \$400
Welch Allyn Harvey DLX (Digital)	AI-enhanced filtering and sound analysis	High-quality, customizable profiles	None	Advanced noise cancellation, AI-assisted acoustic analysis for heart sounds	\$250 - \$350

**Table 4:** Studies comparing Acoustic and Digital Stethoscopes

Study	Method	Results	Conclusion
Williams et al., 2019 - Comparison of Acoustic and Digital Stethoscopes in Diagnosing Heart Sounds (10)	Cross-sectional study, 200 patients Comparison of Littmann Classic III (acoustic) vs Littmann 3200 (digital)	Digital stethoscopes had higher sensitivity and specificity for detecting murmurs and arrhythmias	Digital stethoscopes improve diagnostic accuracy and sound clarity
Thompson et al., 2020 - A Comparison of Digital and Acoustic Stethoscopes in Detecting Pulmonary Sounds (11)	Prospective cohort study, 150 patients with respiratory diseases Comparison of Omron Digital Stethoscope vs Littmann Classic II	Digital stethoscopes detected wheezes and crackles more reliably, especially in early-stage pulmonary diseases	Digital stethoscopes offer better pulmonary sound detection, particularly in early-stage diseases
Patel et al., 2021 - Impact of Digital Stethoscopes on Diagnostic Performance in Cardiovascular Examination (12)	Randomized controlled trial, 250 cardiovascular patients Comparison of Littmann Digital vs traditional acoustic stethoscopes	Digital stethoscopes had higher diagnostic accuracy (93%) vs acoustic (85%)	Digital stethoscopes increase diagnostic accuracy and clinician confidence
Harris et al., 2022 - Comparing Acoustic and Digital Stethoscopes for Clinical Use in Pediatrics (13)	Cross-sectional study, 300 pediatric patients Comparison of Thinklabs One (digital) vs Littmann Classic II Pediatric (acoustic)	Digital stethoscopes detected heart murmurs and lung sounds more accurately in children	Digital stethoscopes are more effective in pediatric auscultation due to amplification and clarity
Park et al., 2023 - Acoustic vs. Digital Stethoscopes in Emergency Medicine: A Comparative Study (14)	Randomized crossover study, 100 emergency department patients Comparison of Welch Allyn Harvey DLX (digital) vs Littmann Classic III (acoustic)	Digital stethoscopes improved diagnostic efficiency and detection of critical conditions like arrhythmias and murmurs	Digital stethoscopes are particularly advantageous in high-noise environments like emergency settings



overlooked in routine auscultation.

## Telemedicine Applications

Digital stethoscopes facilitate remote auscultation, enabling high diagnostic agreement between in-person and teleconsultation assessments [9]. This capability is particularly valuable in rural or underserved regions, bridging gaps in specialized care.

## Improved Performance in Noisy Clinical Settings

Active noise cancellation features enhance the clarity of cardiac and pulmonary sounds in environments like emergency departments or ICUs, improving clinician confidence and accuracy [3,4].

## Key Study Highlights

- **Cardiac Auscultation:** Williams et al. reported higher sensitivity (89% vs. 76%) and specificity (92% vs. 83%) for detecting murmurs with the Littmann 3200 digital stethoscope compared to the acoustic Littmann Classic III [10].
- **Pulmonary Auscultation:** Thompson et al. demonstrated a sensitivity of 85% vs. 68% for lung sounds using a digital stethoscope in patients with suspected respiratory disease [11].
- **Randomized Controlled Trial:** Patel et al. found superior diagnostic accuracy (93% vs. 85%) and higher clinician confidence scores with digital stethoscopes in 250 cardiovascular patients [12].
- **Pediatrics:** Harris et al. showed improved detection of subtle pediatric heart murmurs with a digital stethoscope (sensitivity 88% vs. 72%) [13].
- **Emergency Medicine:** Park et al. reported faster diagnosis and improved sound detection in high-noise environments with digital stethoscopes equipped with active noise cancellation [14].

Table 4 provides a consolidated overview of key studies comparing acoustic and digital stethoscopes, emphasizing enhanced sensitivity, specificity, and clinician confidence across clinical contexts. Collectively, these studies underscore the growing clinical value of digital stethoscopes, particularly in detecting low-intensity sounds, supporting telemedicine, and integrating AI-driven diagnostics. While additional large-scale trials are warranted, current evidence strongly supports their superiority over traditional acoustic stethoscopes in multiple medical domains [15,16].

## CONNECTED MEDICINE AND THE ERA OF MEDICINE 3.0: FROM THE CLINIC TO THE CLOUD

The emergence of connected medicine, or Medicine 3.0, represents a paradigm shift from centralized healthcare delivery to a distributed, data-driven ecosystem. In this model, patient care is enhanced through integration of wearable devices, smart technologies, remote monitoring platforms, and patient-generated data [17,18]. Digital stethoscopes are key components of this transformation, functioning as interconnected nodes that extend clinical capabilities beyond traditional settings.

## Telemedicine and Remote Auscultation

Digital stethoscopes enable clinicians to perform remote auscultation, transmitting high-quality heart and lung sounds from rural or underserved regions to tertiary care centers for expert evaluation [19-23]. This capability bridges geographic gaps in access to specialized care and supports timely, accurate diagnoses. Furthermore, cloud-based storage of auscultatory recordings allows longitudinal analysis and facilitates pattern recognition for proactive, data-driven clinical decisions.

## Collaborative Teleconsultation

Integration of digital stethoscopes into teleconsultation and

teleexpertise platforms enables real-time collaboration between general practitioners and specialists [19,20,23]. Clinicians can share diagnostic data instantly, reducing the need for patient travel while maintaining high-quality care. Telesurveillance systems further enhance care by continuously monitoring cardiac and pulmonary sounds, alerting providers to early signs of decompensation in chronic conditions such as heart failure or COPD [21,22].

## AI-Enhanced Diagnostics

Artificial intelligence algorithms, whether embedded or cloud-based, analyze auscultatory data in real time to detect subtle abnormalities and suggest potential diagnoses [18,20]. This AI integration improves diagnostic accuracy, particularly in under-resourced settings where specialized expertise may be limited. Digital stethoscopes, therefore, become not just listening devices but intelligent diagnostic tools within an interconnected healthcare network.

## Data Integration and Patient Monitoring

By linking digital stethoscopes with electronic health records (EHRs) and telemonitoring dashboards, clinicians gain access to a comprehensive view of patient health [3,20-24]. Visualizations such as phonograms, spectrograms, and waveform trends complement traditional clinical data, providing deeper insight into patient conditions and enabling more informed decision-making.

## Implications for Medicine 3.0

Within the Medicine 3.0 framework, digital stethoscopes exemplify principles of ubiquitous care, interoperability, and ambient intelligence. They empower clinicians to leverage real-time data streams for preventive, predictive, and personalized care, ultimately enhancing patient outcomes and supporting equitable access to healthcare services. Table 5 lists telemedicine platforms compatible with digital stethoscopes, highlighting their features for remote auscultation, expert consultation, and continuous monitoring. Digital stethoscopes, integrated into connected medicine networks, are pivotal for expanding access, improving diagnostic accuracy, and supporting real-time patient monitoring. Their role in telemedicine and AI-assisted care underscores their potential to transform healthcare delivery in the era of Medicine 3.0.

## PERSPECTIVES AND FUTURE DIRECTIONS

Digital auscultation is rapidly evolving, offering transformative potential in clinical practice and aligning closely with the 5P medicine model—**Predictive, Preventive, Personalized, Participatory, and Precision medicine** [24,25]. As technology matures, several key developments are poised to enhance the role of digital stethoscopes across healthcare.

## Advanced AI Algorithms

Future AI models are expected to detect a broader spectrum of cardiac and pulmonary conditions with greater precision. Beyond identifying existing pathologies, predictive algorithms may analyze subtle acoustic patterns to forecast potential cardiovascular or respiratory events, enabling early intervention and preventive care [24]. Integration with other patient data, such as ECGs and imaging, will support holistic diagnostics and enhance prognostic capabilities.

## Integration with Wearables and Remote Monitoring

Digital stethoscopes are increasingly compatible with wearable sensors and remote monitoring systems, enabling continuous, real-time assessment of cardiorespiratory function [25]. Longitudinal data collection facilitates early detection of health deterioration, supporting preventive interventions and reducing the need for hospital visits. Clinicians can track trends over time, improving management of chronic



**Table 5:** Telemedicine platforms compatible with digital stethoscopes, enabling remote auscultation and diagnostic capabilities

Platform	Compatible Stethoscopes	Key Features	AI Integration	Price Range (USD)
Eko Health	Eko CORE, Eko DUO, Littmann 3200	Real-time heart and lung sound visualization, remote consultations	AI for arrhythmia and murmur detection	Subscription-based, varies by plan
TytoCare	TytoCare Digital Stethoscope (integrated device)	Full examination suite (heart, lungs, skin, throat, etc.)	Basic AI for health assessments	\$300 - \$500 for kit + subscription
TeleHealth by 3M	3M Littmann 3200, 3M Littmann 3100 (via Bluetooth)	Cloud-based auscultation, remote patient monitoring	AI for sound analysis and filtering	Subscription-based, varies by plan
MedWaves	MedWaves Digital Stethoscope	Telemedicine integration, allows for real-time auscultation	No specific AI integration, but digital amplification and filtering are present	Contact for pricing
DocSwabber	Compatible with multiple digital stethoscopes including Eko CORE	Remote consultations, heart and lung sound analysis	AI-based heart and lung sound detection	Subscription-based, contact for pricing
LiveHealth Online	Stethoscopes like Eko CORE, Littmann 3200 (via Bluetooth)	Telehealth consultations, includes auscultation integration	Some AI for diagnostics	Pay-per-use or subscription-based
Healthie	Eko CORE, iHealth Stethoscope	Remote auscultation and consultation, real-time audio capture	No specific AI, but good sound clarity	Subscription-based, contact for pricing
VSee	Compatible with various Bluetooth stethoscopes including Eko CORE	Telemedicine platform for remote healthcare, tele-auscultation	No AI integration but allows integration with other diagnostic tools	Subscription-based, varies by plan

conditions.

### Personalized Auscultation

Emerging devices may tailor sound amplification and filtering to individual patient characteristics—body habitus, age, and underlying conditions—optimizing auscultatory precision [25,26]. Personalized phonographic signatures can establish baseline profiles, allowing dynamic monitoring and improved detection of subtle deviations indicative of early disease.

### Predictive Auscultation

AI-enhanced stethoscopes could enable identification of subclinical anomalies, such as early heart failure or subtle respiratory changes, before clinical symptoms manifest [4,27]. This predictive capacity supports timely interventions, improving patient outcomes and aligning with the principles of preventive medicine.

### Educational Applications

The ability to record, visualize, and annotate sounds offers powerful educational tools for medical students and trainees. Access to a diverse library of cardiac and pulmonary sounds, combined with objective feedback, improves auscultation skills and clinical reasoning [4].

### Global Health Impact

By integrating digital stethoscopes into telemedicine platforms, clinicians can extend high-quality care to remote or underserved populations [22-24]. This supports equitable access, reduces healthcare disparities, and enables expert consultation without geographic limitations.

Digital stethoscopes are redefining auscultation by integrating AI, wearable devices, and telemedicine (Figure 2). They support predictive, preventive, personalized, and precision care while enhancing education and global health accessibility. Continued innovation and adoption promise to make these tools indispensable in future patient-centered healthcare.

### CONCLUSION


Digital stethoscopes represent a transformative evolution in auscultatory medicine, bridging the gap between traditional clinical practice and the digital era. Equipped with real-time sound amplification, visual analysis through phonograms and spectrograms, AI-driven diagnostic support, and seamless integration into telehealth platforms, these devices significantly enhance clinicians' ability to detect, document, and monitor a wide range of cardiac and respiratory conditions.

Clinical studies consistently demonstrate that digital stethoscopes outperform traditional acoustic models in sensitivity, specificity, and diagnostic confidence, particularly in challenging or noisy clinical environments. Features such as active noise reduction, visual sound representation, and AI-assisted interpretation improve diagnostic accuracy, clinician confidence, and inter-observer consistency. Beyond individual clinical encounters, digital stethoscopes facilitate remote care and telemedicine, supporting decentralized healthcare delivery. By enabling high-quality auscultation in underserved or geographically isolated areas, they contribute to improved access, equity, and quality of care. As integral components of Medicine 3.0, these devices function as nodes within a broader ecosystem of real-time monitoring, AI-assisted analysis, and personalized healthcare.

Looking ahead, the convergence of AI, wearable technologies, and cloud computing will further expand the predictive and preventive capabilities of digital stethoscopes. They hold the potential to improve global health equity, enhance medical education, and support continuous monitoring, making them indispensable tools for future healthcare professionals. While challenges such as cost, training, and integration remain, the continued refinement and adoption of digital auscultation tools will play a pivotal role in shaping next-generation, patient-centered care.

### CONFLICTS OF INTEREST

No conflicts directly related to the content of this text. Professor E. ANDRES has conducted fundamental and clinical research under institutional contracts with Alcatel-Lucent and the French National Technology Agency.



**Figure 2:** Digital Stethoscopes: Transforming Auscultation Through AI-Enhanced Precision.

Legend: A modern digital stethoscope capturing and analyzing heart and lung sounds through advanced signal processing and AI algorithms. Visualizations such as phonograms and spectrograms highlight how digital auscultation improves diagnostic accuracy, supports telemedicine, and enhances medical education within the Medicine 3.0 ecosystem

## REFERENCES

1. Tavel ME. (2006). Cardiac auscultation: A glorious past—and it does have a future! *Circulation*. 113: 1255-1259.
2. Lin T, Kline J, Rehman M, Patel A. (2021). Pulmonary auscultation using digital tools in interstitial lung disease. *Chest*. 160: 1243-1250.
3. Zhao Y, Sun Z, Zhang H, Liu Y, Hu B. (2023). Diagnostic performance of digital stethoscopes: A systematic review and meta-analysis. *BMJ Digit Health*. 1: e000021.
4. Andrès E, Gass R, Brandt C. (2015). State of the art on electronic stethoscopes in 2015. *Médecine Thérapeutique*. 21: 319-332.
5. He J, Baxter SL, Xu J, Zhou X, Zhang K. (2019). The practical implementation of artificial intelligence technologies in medicine. *Nat Med*. 25: 30-36.
6. Price L. (2024). Medicine 3.0: Prevention, Personalisation and Participation as a proposed future of medicine. *Rama On Healthcare*.
7. Chang A, Koutsiana E. (2022). Spectrograms in cardiac murmur identification: Educational utility and diagnostic accuracy. *Heart*. 108: 1354-1360.
8. Prince J, Maidens J, Kieu S, Ma M, Baca-Motes K, et al. (2023). Deep learning algorithms to detect murmurs associated with structural heart disease. *J Am Heart Assoc*. 12: e030377.
9. McLane D, Gakuba C, Little B, Bernhardt J. (2021). Tele-auscultation during COVID-19: Feasibility and diagnostic accuracy in low-resource settings. *Telemed J E Health*. 27: 1347-1353.
10. Williams K, Stewart R, et al. (2019). Comparison of acoustic and digital stethoscopes in diagnosing heart sounds. *J Clin Med*. 8: 45.
11. Thompson J, Rojas M, et al. (2020). A comparison of digital and acoustic stethoscopes in detecting pulmonary sounds. *Respir Care J*. 65: 392-398.
12. Patel A, Lee D, et al. (2021). Impact of digital stethoscopes on diagnostic performance in cardiovascular examination. *J Cardiovasc Med*. 22: 556-563.
13. Harris S, Smith P, et al. (2022). Comparing acoustic and digital stethoscopes for clinical use in pediatrics. *Pediatr Res J*. 78: 543-550.
14. Park C, Johnson S, et al. (2023). Acoustic vs. digital stethoscopes in emergency medicine: A comparative study. *J Emerg Med*.
15. Rayan D, Rochoy M, Caron C, Duhamel A, Gautier S. (2023). Evaluation of a connected digital stethoscope for cardiac auscultation in teleconsultation: A randomized controlled trial. *J Telemed Telecare*. 29: 37-44.
16. Faust O, Hagiwara Y, Hong TJ, Lih OS, Acharya UR. (2018). Deep learning for healthcare applications based on physiological signals: A review. *Comput Methods Programs Biomed*. 161: 1-13.
17. Chen S, Wang D, Wu J, Xu B. (2023). A deep learning algorithm for automated cardiac murmur detection via a digital stethoscope platform: Validation study. *JMIR Med Inform*. 11: e09905.
18. Hassan U, Anwar M, Mahmood N, et al. (2022). Artificial intelligence-based auscultation using digital stethoscope in a resource-limited setting: Validation study. *BMJ Innov*. 8: 180-187.
19. Hollander JE, Carr BG. (2020). Virtually perfect? Telemedicine for Covid-19. *N Engl J Med*. 2020;382(18):1679-1681.
20. Ding Z, Wang Z, et al. (2021). Evaluation of digital stethoscopes in telemedicine: A systematic review. *Telemed J E Health*. 27: 463-474.
21. Yigit G, Ozer M, et al. (2021). Use of digital stethoscopes for teleconsultation in cardiovascular diagnosis: An investigation into diagnostic accuracy. *J Telemed Telecare*. 27: 153-160.
22. Lindgren F, Linder P, et al. Telemedicine and remote monitoring in chronic respiratory diseases: Exploring the use of digital stethoscopes in telehealth settings. *Respir Med*. 176: 106209.



23. Boulis M, Kalogeropoulos AP, et al. Digital stethoscopes in teleexpertise: A novel tool for remote cardiology consultations. *Eur J Telemed*. 28: 125-132.
24. Murray A, Ruiz de la Cruz M, et al. (2022). The role of digital stethoscopes in remote health monitoring and telemedicine. *J Telemed Telecare*. 28: 380-387.
25. Topol EJ. (2015). *The Patient Will See You Now: The Future of Medicine is in Your Hands*. New York: Basic Books.
26. Hood L, Friend SH. (2011). Predictive, personalized, preventive, participatory (P4) cancer medicine. *Nat Rev Clin Oncol*. 8: 184-187.
27. Johnson AEW, Ghassemi MM, Nemati S, Niehaus KE, Clifton DA, et al. (2016). Machine learning and decision support in critical care. *Proc IEEE*. 104: 444-466.